

HAMMERTECH

Reaming

Implant Kit

UP TO **200%** more bone to bone contact

Uniform articular tissue removal

Infinite degrees of freedom



Hemispherical reaming design provides more bone to bone contact surface area compared to traditional linear cuts

Matching ball and socket geometry provides infinite degrees of freedom while maintaining joint apposition

Designed to maintain toe length compared to linear cuts



Ordering Information

Description	Diameter	Part Number
8mm Reaming Implant Kit	8mm	HT-60-0808-S
10mm Reaming Implant Kit	10mm	HT-60-1010-S
12mm Reaming Implant Kit	12mm	HT-60-1212-S

FUSION ORTHOPEDICS, LLC.
4135 S. Power Rd. Suite 110
Mesa, Arizona 85212 USA
(800)-403-6876
Sales@fusionorthopedics.com

Fusion
ORTHOPEDICS

FDA Cleared.
All Rights Reserved. HammerTech is a
Trademark of Fusion Orthopedics, LLC.
D-4060-SLSX-2 Rev A

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Fusion Orthopedics does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Fusion's product offerings. A surgeon must always refer to the package insert, product label and / or instructions for use before using any Fusion Orthopedics product.