

TALARLIFT

STS

Subtalar System



TALARLIFT

STS

Subtalar System



Cylindrical



Conical



Anatomical

We Believe in Everything Better

- **3 Implant Styles In One Tray** and 18 unique implants make TalarLift the most comprehensive subtalar offering on the market.
- **Progressively Softened Threads** are designed to minimize the edge effects that may lead to pain.
- **Apertures** allow for soft tissue ingrowth to resist lateral migration of the implant.
- **Multiple Size Options** are designed for varying patient anatomy with length and diameters that scale with size.
- **Supporting Instrumentation** features color-coded trials for verification of placement and ease of use.

IMPLANT ORDERING



	Cylindrical	Conical	Anatomical
7mm	TL-CY-0007	TL-CN-0007	TL-AT-0007
8mm	TL-CY-0008	TL-CN-0008	TL-AT-0008
9mm	TL-CY-0009	TL-CN-0009	TL-AT-0009
10mm	TL-CY-0010	TL-CN-0010	TL-AT-0010
11mm	TL-CY-0011	TL-CN-0011	TL-AT-0011

FUSION ORTHOPEDICS, LLC.

4135 S. Power Rd. Suite 110
Mesa, Arizona 85212 USA
(800)-403-6876
Sales@fusionorthopedics.com



FDA Cleared

All Rights Reserved. TalarLift is
a Trademark of Fusion Orthopedics, LLC.
D-4560-SLSX Rev A

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Fusion Orthopedics does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Fusion's product offerings. A surgeon must always refer to the package insert, product label and / or instructions for use before using any Fusion Orthopedics product.